## EN American Eagle Hand Instruments Reprocessing Instructions

### AEI Dental Hand Instruments and Accessories

New AEI Dental Hand Instruments must be sterilized prior to patient use. The instruments are reusable and should be sterilized using steam sterilization only.

### Product Description

The AEI Dental Hand Instrument product Line includes instrument groups of diagnostic hand instruments, periodontal hand instruments, and restorative hand instruments.

### Intended Use

AEI Dental Hand Instruments are only to be used in professional dental settings, by trained dental professionals and are to be used in accordance with industry standard dental practices for the instrument’s intended purpose. AEI XP Dental Hand Instruments are “Sharpen Free” and are not intended to be sharpened. The XP cutting edge protection and non-stick abilities will be negated if sharpened or damaged.

### Warnings and Precautions

- Improper use of AEI Dental Hand Instruments is unsafe and may cause harm to the dental professional and/ or patient.
- AEI Dental Hand Instruments are resistant to reprocessing and end of life is determined by wear and damage from use.

- Prior to every use of AEI Dental Hand Instruments, the clinician must: (1) Clean off visible debris from the instrument. (2) Place in a decontamination/disinfection bath per the attached Decontamination and Disinfection procedures listed below. (3) Use steam sterilization via an autoclave. Additional instructions on each of these steps can be found below. Prior to any additional use, the clinician must first visually inspect the AEI Dental Hand Instrument(s) for any damage, corrosion, or defects that could cause the device to be unsuitable for patient use. If any of these characteristics are found, discard the device in accordance with local regulations for possible contaminated devices if not sterilized.

- Use only FDA approved medical or dental reprocessing equipment and neutral PH enzymatic cleaners. The provided reprocessing instructions have been validated. It is the user’s responsibility to follow the manufacturer’s instructions for alternative reprocessing procedures.

- Use of cleaning chemicals that contain Aqua Regia, iron chloride, sulfuric and hydrochloric acid, and iodine should be avoided. See the AEI website for a full list of compatible and incompatible cleaners and chemicals.

- It is the users responsibility to verify the functionality of a disinfection unit and steam autoclaves at regular intervals by having it serviced and calibrated according to the manufacturer’s recommendations.

- If a device arrives unsuitable for patient use, contact AEI Customer Service. Small tears and cuts to the instrument packaging is acceptable if the instrument is undamaged.

### Reprocessing Instructions

- Disassemble all devices prior to sterilization. Quik-Tips should be removed from the handle and separated for proper sterilization.

- Disassemble all devices prior to manual cleaning and disinfection. Quik-Tips should be removed from the handle and separated for proper manual cleaning and disinfection.

- The automatic washing/disinfection reprocessing efficacy was validated with the Quik-Tips assembled. Only the automatic washing/disinfection method does not require the Quik-tips to be separated from their handles.

- Do not use an abrasive pad, metal brush, or steel wool for cleaning as they could cause micro-abrasions that could damage the AEI Dental Hand Instrument.

- Dissimilar metals should not be processed together. If unlike metals are batched together, a potential for electrolytic reactions between dissimilar metals will exist, which can produce etching and corrosion on the instrument surfaces.

- The devices should not contact each other to prevent contact damage and surface shielding of cleaners or steam.

### Point of Use Debris Removal (Optional)

- Remove bioburden and debris from devices as soon as possible to prevent residual adhesion.

- If manual decontamination cannot be performed immediately, place the device

- in a covered container, submerged in distilled water or in a bath of neutral PH enzymic cleaner, follow the cleaner manufactures instructions. Do not soak longer than 60 minutes.

### Manual Decontamination

- Using a sink dedicated to contaminated device handling, fully dissemble the device and rinse in cold tap water for 30 seconds to remove excess visible debris from the instrument.

- In a tub, prepare a fresh bath of cleaning solution with a PH neutral enzymatic cleaner, follow the manufacturer instructions for concentration, time, and water quality.

- Submerge the instruments and use a soft-bristled nylon brush on each device for at least 30 seconds hitting all surfaces of the device. See warnings in Reprocessing Instructions section for authorized brush materials.

- Thoroughly rinse the instrument in tap water for 30 seconds to remove visible contaminants and cleaning solution.

- Blot the device dry with a clean, non-shedding wipe, inspecting the device to ensure all contaminants have been removed. Repeat as necessary if visible debris or contamination is still visible.

### Manual Disinfection

- Prepare a fresh sonic bath with a cleaner/disinfectant, follow the manufacturer instructions for concentration, time, and water quality.

- Load instruments into a cassette or ultrasonic mesh basket. See warnings in Reprocessing Instructions section for instrument contact and reprocessing of dissimilar metals.

- Process the devices, following the manufacturer instructions for concentration, time, temperature, and water specification. Sonicate for 15 minutes if no specifications given.

- Remove the device from the sonic bath, for best results rinse in distilled or filtered water for 30 seconds to remove any cleaning fluids.

- Allow the device to dry following the rinse, using a clean non-shedding wipe if necessary to expedite drying.

- Place the disinfected instruments in a disinfected cassette or sterilization pouch for sterilization.

### Auto Washer Cleaning/Disinfection Settings

- Always follow the Auto Washer/Disinfectant and sterilizing equipment manufacturer’s recommendations for time, temperature, and water quality.

- Carefully load the instruments into a cassette or other container. See warnings in Reprocessing Instructions section for instrument contact and reprocessing of dissimilar metals.

- AEI recommends the following automatic washer-disinfectior cycle: Pre-rinse with cold tap water for 1 minute, Wash at 55°C for 5 minutes with distilled water and 0.5% Neodisher Mediclean Forte, Post rinse in cold demineralized water for 2 minutes, Thermal Disinfection at 93°C for 5 minutes, use machine settings for drying.

- Place disinfected instruments into a sterilizing pouch, cassette, or other disinfected container for transfer to sterilization.

### Sterilizing Settings

- Place instruments into the Dynamic Air Removal steam autoclave unit and process the instruments according to the manufacturer’s recommendations, but minimally in accordance with AEI’s validated Dynamic Air Removal instrument sterilization procedure below for wrapped instruments. A sterilization Test Strip is recommended to be placed in each autoclave cycle to monitor the sterilization process. The Test Strip will not record the need for routine maintenance and sterilization efficacy checks.

- See Reprocessing Instructions section for warnings of instrument contact and reprocessing of dissimilar metal.

	Country	Temperature	Time	Drying Time
	Europe	134°C (273°F)	≥ 3 Min.	30 Min.
	USA	132°C (270°F)	≥ 4 Min.	30 Min.

- Remove the instruments from the autoclave as soon as possible following completion of the sterilization-drying cycle.

- After sterilization, confirm that the instruments are fully dry. Store the instruments in appropriate packaging that indicates the status of the sterilization and maintains the instrument sterilization until used. Do not store with non-sterilized instruments.

*Thank you for trusting AEI for your dental hand instrument needs.*

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