



Compatible Solutions List

These solutions have shown to be compatible with AEI instruments. Consult AEI IFU for reprocessing instructions to follow cleaning and sterilizing times and temperatures.

PRODUCT	MANUFACTURER	TEST RESULTS
Monarch Enzymatic Solution (without SXS)	Air Techniques	Compatible
OptiSonic™ Super-Concentrated (General Purpose Ultrasonic Cleaning Solution)	Benco	Compatible
Benco Enzyme Tablets	Benco	Compatible
Benco General Purpose Ultrasonic Cleaning Solution	Benco	Compatible
BiozymeLT Enzyme Detergent	Biotrol	Compatible
BioSonic Enzymatic Ultrasonic Cleaner	Coltene/Whaledent Inc	Compatible
BrandMax Triple Enzyme Cleaner	BrandMax	Compatible
Cetylite® Power Cleaner	Cetylite Industries Inc	Compatible
Crosstex® General Purpose Ultrasonic Cleaner	Cetylite Industries Inc	Compatible
Cetyl-Zyme Pro-Am Dual Enzymatic Detergent Concentrate	Cetylite Industries Inc	Compatible
Cetyl-Zyme Pro-Am Foam	Crosstex	Compatible
Instrument Forte	Helvemed	Compatible
Enzymax® Liquid	Hu-Friedy	Compatible
Enzymax® Powder	Hu-Friedy	Compatible
Enzymax Earth® Tablets	Hu-Friedy	Compatible
Patterson General Purpose Type1	Patterson	Compatible
Pearson Quality Super	Pearson	Compatible
Brite Sheild	Pemier	Compatible
ReSURGE™ Instrument Cleaning Solution	Sultan Healthcare	Compatible
Clean & Simple Tablets	Tuttnauer	Compatible
GZYME3	Germiphene Corp.	Compatible
Vigilance	Harry J. Bosworth Co.	Compatible
Citrizyme	Pascal International, Inc.	Compatible
Good Vibrations	Certol International LLC	Compatible
Resurge	Sultan Healthcare	Compatible
Empower Dual Enzymatic Detergent	Metrex Corporation	Compatible
Eco 4 Neutral All-Purpose Cleaner	Sonix4	Compatible
Purit Cide It	Biotrol	Compatible
Gleam II Ultrasonic Cleaning Solution	Germiphene	Compatible
Microbex Microvac Cleaning Solution	Germiphene	Compatible
Supersonic All Purpose Cleaning Solution	Germiphene	Compatible
Z3™ Multi-Purpose Enzymatic Detergent	Benco	Incompatible
Opti•Cide3	Biotrol	Incompatible
neodisher® FA	Chemische Fabrik Dr. Weigert GmbH&Co. KG	Incompatible
neodisher® MA	Chemische Fabrik Dr. Weigert GmbH&Co. KG	Incompatible
neodisher® Mielclear	Chemische Fabrik Dr. Weigert GmbH&Co. KG	Incompatible
neodisher® N	Chemische Fabrik Dr. Weigert GmbH&Co. KG	Incompatible
Denta Zyme®	Ecolab	Incompatible
COEZYME™ Enzymatic Detergent	GC America	Incompatible
GZYMXE	Germiphene Corp	Incompatible
VibraKleen® E Tablets	Heraeus Kulzer	Incompatible
EZ-Zyme® All Purpose Cleaner	Integra Miltex	Incompatible
Multi-Purpose Enzyme	Patterson	Incompatible
Ultra Dose General Purpose Cleaner	L&R Manufacturing Company	Incompatible
Microbex Microzyme	Germiphene Corp.	Incompatible